

Linezolid

Catalog No: #S0391


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Linezolid
Brief Description	Inhibitors
Purification	99.00%
Target Name	MAO inhibitor
Calculated MW	337.35
Formulation	C16H20FN3O4
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:ImmunologySMILES:CC(=O)NC[C@H]1CN(C(=O)O1)C2=CC(=C(C=C2)N3CCOCC3)FPathways:MetabolismReceptor:AOA;
AOBBoling pt:585.5C?Melting pt:176-178CSolubility:Ethanol: 8 mg/mL (23.71 mM); DMSO: 67 mg/mL (198.6 mM)Appearance:White
PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be
stored below -20 °C for several month.

References

1. Stevens DL, et al. Expert Rev Anti Infect Ther. 2004 Feb;2(1):51-9.

Note: This product is for in vitro research use only