

## Moxisylyte hydrochloride

Catalog No: #S0394


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Moxisylyte hydrochloride
Brief Description	Inhibitors
Purification	99.00%
Target Name	Adrenergic Receptor inhibitor
Calculated MW	315.84
Formulation	C16H26ClNO3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Cardiovascular System  
 SMILES: Cl.CC(C)C1=C(OCCN(C)C)C=C(C)C(OC(C)=O)=C1  
 Pathways: Neuroscience; GPCR/G  
 Protein Receptor: Adrenergic Receptor  
 Boiling pt: 371°C  
 Melting pt: Solubility: DMSO: 10 mM  
 Appearance: White Powder  
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Imagawa A, et al. Life Sci. 1989;44(9):619-23.

Note: This product is for in vitro research use only