

# Pyrantel tartrate

Catalog No: #S0395


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Pyrantel tartrate
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	356.39
Formulation	C15H20N2O6S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Other  
 SMILES: CN1CCCN=C1/C=C/C2=CC=CS2.[C@@H]([C@H](C(=O)O)O)(C(=O)O)O  
 Pathways: Others  
 Receptor: Others  
 Boiling pt: 324.4°C  
 Melting pt: 148-150°C  
 Solubility: DMSO: 10 mM  
 Appearance: Off-white Color Powder  
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Pratt SE, et al. Am J Vet Res. 1981 May;42(5):871-2.

Note: This product is for in vitro research use only