

Pipemidic acid

Catalog No: #S0396


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Pipemidic acid
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	303.32
Formulation	C ₁₄ H ₁₇ N ₅ O ₃
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Immunology
 SMILES: CCN1C=C(C(O)=O)C(=O)C2=CN=C(N=C12)N1CCNCC1
 Pathways: Others
 Receptor: Others
 Boiling pt: Melting pt: 251-255°C
 Solubility: DMSO: 10 mM
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Hoeglund AB, et al.
 J Med Chem. 2010 Feb 11;53(3):1056-66.

Note: This product is for in vitro research use only