

Carzenide

Catalog No: #S0398


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Carzenide
Brief Description	Inhibitors
Purification	98.65%
Target Name	Others
Calculated MW	201.2
Formulation	C7H7NO4S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Nervous system
 SMILES: NS(=O)(=O)C1=CC=C(C=C1)C(=O)O
 Pathways: Others
 Receptor: Others
 Boiling pt: Melting pt: 285-295°C
 Solubility: DMSO: 10 mM
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Petzinger E. Dtsch Tierarztl Wochenschr. 1992 Oct;99(10):412-5.

Note: This product is for in vitro research use only