

## Methyl salicylate

Catalog No: #S0399


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Methyl salicylate
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	152.15
Formulation	C <sub>8</sub> H <sub>8</sub> O <sub>3</sub>
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Immunology  
SMILES: COC(=O)C1=CC=CC=C1O  
Pathways: Others  
Receptor: Others  
Boiling pt: 222°C  
Melting pt: -8°C  
Solubility: Soluble in Ethanol, slightly soluble in Water  
Appearance: Liquid  
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Spiandore M, et al. Drug Test Anal. 2014 Jun;6 Suppl 1:67-73.

Note: This product is for in vitro research use only