

## Octocrilene

Catalog No: #S0402


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Octocrilene
Brief Description	Inhibitors
Purification	98.26%
Target Name	Others
Calculated MW	361.48
Formulation	C24H27NO2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Other  
 SMILES: CCCCC(CC)COC(=O)C(C#N)=C(C1=CC=CC=C1)C1=CC=CC=C1  
 Pathways: Others  
 Receptor: Others  
 Boiling pt: 218°C  
 Melting pt: -10°C  
 Solubility: DMSO: 10 mM  
 Appearance: Liquid  
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

- Hanson KM. Free Radic Biol Med. 2006 Oct 15;41(8):1205-12.

Note: This product is for in vitro research use only