

Hydroxyquinoline

Catalog No: #S0404


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Hydroxyquinoline
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others inhibitor
Calculated MW	145.16
Formulation	C ₉ H ₇ NO
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Immunology
SMILES: OC1=CC=CC2=C1N=CC=C2
Pathways: Others
Receptor: Others
Boiling pt: 267°C
Melting pt: 70-73°C
Solubility: Water: 556 mg/l at 20 °C.
Appearance: Gray Crystal
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.
Stock solution can be stored below -20 °C for several month.

References

1. Freitas LB, et al. Eur J Med Chem. 2014 Sep 12;84:595-604.

Note: This product is for in vitro research use only