

Dehydrocholic acid

Catalog No: #S0408


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Dehydrocholic acid
Brief Description	Inhibitors
Purification	95.00%
Target Name	Others
Calculated MW	402.52
Formulation	C ₂₄ H ₃₄ O ₅
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Metabolism

system SMILES: C[C@H](CCC(=O)O)[C@H]1CC[C@@H]2[C@@]1(C(=O)C[C@H]3[C@H]2C(=O)C[C@H]4[C@@]3(CCC(=O)C4)C)C
 Pathways: Others
 Receptor: Others
 Boiling pt: 581.5°C
 Melting pt: 238-240°C
 Solubility: DMSO: 10 mM
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Mohamed AA, Matijevi? E. J Colloid Interface Sci. 2012 Feb 15;368(1):625-8.

Note: This product is for in vitro research use only