

Phthalylsulfacetamide


Catalog No: #S0409

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Phthalylsulfacetamide
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	362.36
Formulation	C16H14N2O6S
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:ImmunologySMILES:CC(=O)NS(=O)(=O)C1=CC=C(NC(=O)C2=CC=CC=C2C(O)=O)C=C1Pathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. ASHER LM. Gastroenterology. 1958 Apr; 34(4):659-60.

Note: This product is for in vitro research use only