

## Calcium gluconate

Catalog No: #S0411


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Calcium gluconate
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	430.37
Formulation	C <sub>12</sub> H <sub>22</sub> CaO <sub>14</sub>
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Metabolism

system SMILES: C([C@H]([C@H]([C@@H]([C@H](C(=O)[O-])O)O)O)O)O.C([C@H]([C@H]([C@@H]([C@H](C(=O)[O-])O)O)O)O).[Ca+2] Pathways: Others Receptor: Others Boiling pt: 673.6 Melting pt: 195 Solubility: Water: 30 mg/mL Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Drop LJ, Cullen DJ. Br J Anaesth. 1980 May;52(5):501-5.

Note: This product is for in vitro research use only