

Cinromide

Catalog No: #S0413


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Cinromide
Brief Description	Inhibitors
Purification	97.73%
Target Name	Others
Calculated MW	254.12
Formulation	C ₁₁ H ₁₂ BrNO
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
 SMILES: CCNC(=O)/C=C/C1=CC(=CC=C1)Br
 Pathways: Others
 Receptor: Others
 Boiling pt: 417.5°C
 Melting pt: 89-91°C
 Solubility: DMSO: 10 mM
 Appearance: White Crystal
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Fromm GH, et al. *Epilepsia*. 1983 Aug; 24(4):394-400.

Note: This product is for in vitro research use only