

Pentetic Acid

Catalog No: #S0420


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Pentetic Acid
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	393.35
Formulation	C ₁₄ H ₂₃ N ₃ O ₁₀
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Metabolism system
 SMILES: OC(=O)CN(CCN(CC(O)=O)CC(O)=O)CCN(CC(O)=O)CC(O)=O
 Pathways: Others
 Receptor: Heavy metal ion
 Boiling pt: Melting pt: 219-220°C
 Solubility: Soluble in Water
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Gi M, et al. Antimicrob Agents ChemOthers. 2014 Dec; 58(12):7205-14.

Note: This product is for in vitro research use only