

Inositol


Catalog No: #S0421

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Inositol
Brief Description	Inhibitors
Purification	99.00%
Target Name	Vitamin antagonist
Calculated MW	180.16
Formulation	C ₆ H ₁₂ O ₆
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Metabolism system
SMILES: OC1C(O)C(O)C(O)C(O)C1O
Pathways: Metabolism
Receptor: VBB
Boiling pt: 291.3°C
Melting pt: 222-227°C
Solubility: Soluble in Water
Appearance: White Crystal
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Chen X, et al. Nucleic Acids Res. 2002 Jan 1; 30(1):412-5.

Note: This product is for in vitro research use only