

## Talc

Catalog No: #S0422


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Talc
Brief Description	Inhibitors
Purification	0.0450mm
Target Name	Others
Calculated MW	379.27
Formulation	3MgO4SiO2H2O
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:ImmunologySMILES:O.O=[Mg].O=[Mg].O=[Mg].O=[Si]=O.O=[Si]=O.O=[Si]=O.O=[Si]=OPathways:OthersReceptor:OthersBoiling pt:Melting pt:800CSolubility:DMSO: 0.1 mg/mL (0.26 mM)Appearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Ahmed L, et al. Chest. 2014 Dec; 146(6):e190-4.

Note: This product is for in vitro research use only