

Edetate?Trisodium


Catalog No: #S0423

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Edetate?Trisodium
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	358.19
Formulation	C ₁₀ H ₁₃ N ₂ Na ₃ O ₈
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Metabolism

systemSMILES:O.O.[Na+].[Na+].[Na+].OC(=O)CN(CCN(CC([O-])=O)CC([O-])=O)CC([O-])=OPathways:OthersReceptor:metal ionBoling

pt:614.2CMelting pt:237CSolubility:Soluble in WaterAppearance:White CrystalRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Stowers JM, et al. Lancet. 1967 Jan 21;1(7482):124-7.

Note: This product is for in vitro research use only