

Bismuth subsalicylat


Catalog No: #S0424

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Bismuth subsalicylat
Brief Description	Inhibitors
Purification	99.00%
Target Name	COX inhibitor
Calculated MW	362.09
Formulation	C7H5BiO4
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Immunology
SMILES: O.O=C1O[Bi]OC2=CC=CC=C12
Pathways: Immunology/Inflammation; Neuroscience
Receptor: COX-1; COX-2
Boiling pt: Melting pt: >350°C
Solubility: DMSO: 0.04 mg/mL (0.11 mM)
Appearance: White Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Pitz AM, et al. Gut Microbes. 2015 Mar 4; 6(2):93-100.

Note: This product is for in vitro research use only