

Bismuth subsalicylat

Catalog No: #S0424


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Bismuth subsalicylat
Brief Description	Inhibitors
Purification	99.00%
Target Name	COX inhibitor
Calculated MW	362.09
Formulation	C7H5BiO4
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:ImmunologySMILES:O.O=C1O[Bi]OC2=CC=CC=C12Pathways:Immunology/Inflammation; NeuroscienceReceptor:COX-1; COX-2Boiling pt:Melting pt:>350CSolubility:DMSO: 0.04 mg/mL (0.11 mM)Appearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Pitz AM, et al. Gut Microbes. 2015 Mar 4; 6(2):93-100.

Note: This product is for in vitro research use only