

COP  $\zeta$ 1 Polyclonal Antibody

Catalog No: #40779

Package Size: #40779-1 50ul #40779-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	COP $\zeta$ 1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms Mk
Specificity	COP $\zeta$ 1 Polyclonal Antibody detects endogenous levels of COP $\zeta$ 1 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the N-terminal region of human COP $\zeta$ 1.
Target Name	COP $\zeta$ 1
Other Names	COPZ1; COPZ; CGI-120; HSPC181; Coatomer subunit zeta-1; Zeta-1-coat protein; Zeta-1 COP
Accession No.	Swiss-Prot: P61923NCBI Gene ID: 22818
Uniprot	P61923
GeneID	22818;
SDS-PAGE MW	20kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/10000.

Not yet tested in other applications.

## Images

Western Blot analysis of COS7 cells using COP  $\zeta$ 1 Polyclonal Antibody

---

Note: This product is for in vitro research use only