

## Diclofenac diethylamine

Catalog No: #S0432


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Diclofenac diethylamine
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others inhibitor
Calculated MW	369.29
Formulation	C18H22Cl2N2O2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Immunology  
 SMILES: CCNCC.O=C(O)CC1=CC=CC=C1NC1=C(Cl)C=CC=C1Cl  
 Pathways: Others  
 Receptor: Others  
 Boiling pt: Melting pt: 145-148°C  
 Solubility: Ethanol: 74 mg/mL (200.38 mM); DMSO: 74 mg/mL (200.38 mM)  
 Appearance: White Powder  
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Dastidar SG, et al. Int J Antimicrob Agents, 2000, 14(3), 249-251.

Note: This product is for in vitro research use only