

## Dicurone

Catalog No: #S0434


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Dicurone
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	176.12
Formulation	C6H8O6
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:Metabolism systemSMILES:C(=O)[C@@H]([C@@H]1[C@@H]([C@@H](C(=O)O1)O)O)OPathways:OthersReceptor:OthersBoiling pt:403.5CMelting pt:172-175CSolubility:DMSO: 10 mMAppearance:White CrystalRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Suzuki S, et al. J Chromatogr Sci. 1998 Jul;36(7):357-60.

Note: This product is for in vitro research use only