

## Sulfinpyrazone

Catalog No: #S0435


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Sulfinpyrazone
Brief Description	Inhibitors
Purification	98.09%
Target Name	COX inhibitor; Others
Calculated MW	404.48
Formulation	C23H20N2O3S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Cardiovascular

System SMILES: O=C1C(CCS(=O)C2=CC=CC=C2)C(=O)N(N1C1=CC=CC=C1)C1=CC=CC=C1 Pathways: Immunology/Inflammation; Others;

Neuroscience Receptor: COX; cMOAT Boling pt: Melting pt: Solubility: Soluble in DMSO Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Pay GF, et al. Haemostasis. 1981;10(3):165-75.

Note: This product is for in vitro research use only