

## Acexamic Acid


Catalog No: #S0436

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Acexamic Acid
Brief Description	Inhibitors
Purification	98.00%
Target Name	Amino Acids and Derivatives inhibitor; COX inhibitor
Calculated MW	173.21
Formulation	C8H15NO3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Other  
SMILES: CC(=O)NCCCCCC(=O)O  
Pathways: Immunology/Inflammation; Others; Neuroscience  
Receptor: Amino Acids and Derivatives; COX  
Boiling pt: Melting pt: 102-104°C  
Solubility: soluble in Water  
Appearance: White Powder  
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Lee U, et al. *Xenobiotica*. 2012 Mar;42(3):310-5.

Note: This product is for in vitro research use only