

L-Tryptophan

Catalog No: #S0439


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	L-Tryptophan
Brief Description	Inhibitors
Purification	98.00%
Target Name	Amino Acids and Derivatives
Calculated MW	204.23
Formulation	C ₁₁ H ₁₂ N ₂ O ₂
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Nervous system
SMILES: C1=CC=C2C(=C1)C(=CN2)C[C@@H](C(=O)O)N
Pathways: Others
Receptor: Amino Acids and Derivatives
Boiling pt: Melting pt: 289-290°C
Solubility: soluble in Water
Appearance: White Crystal
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Tsuchiya W, et al. Nucleic Acids Symp Ser (Oxf). 2004;(48):185-6.

Note: This product is for in vitro research use only