

## Pasiniiazid

Catalog No: #S0441


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Pasiniiazid
Brief Description	Inhibitors
Purification	97.00%
Target Name	COX; NADPH
Calculated MW	290.28
Formulation	C13H14N4O4
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Immunology  
 SMILES: NNC(=O)C1=CC=NC=C1.NC1=CC(O)=C(C=C1)C(O)=O  
 Pathways: Immunology/Inflammation; Metabolism; Neuroscience  
 Receptor: COX; NADPH  
 Boiling pt: 380.8°C  
 Melting pt: Solubility: Soluble in Water  
 Appearance: Brownness Powder  
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Pym AS, et al. Mol Microbiol. 2001 May;40(4):879-89.

Note: This product is for in vitro research use only