

Cinchophen


Catalog No: #S0456

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cinchophen
Brief Description	Inhibitors
Purification	99.00%
Target Name	COX inhibitor
Calculated MW	249.26
Formulation	C16H11NO2
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Nervous systemSMILES:OC(=O)C1=CC(=NC2=CC=CC=C2)C1=CC=CC=C1Pathways:Immunology/Inflammation;
NeuroscienceReceptor:COXBoling pt:456.9CMelting pt:214-215CSolubility:Ethanol: 5 mg/mL (20.05 mM); DMSO: 50 mg/mL (200.59
mM)Appearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a
while.Stock solution can be stored below -20 °C for several month.

References

1. Cutrn Prieto C, et al. Med Clin (Barc). 1991 Jun 15; 97(3):104-6.

Note: This product is for in vitro research use only