

Ellagic acid

Catalog No: #S0465


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Ellagic acid
Brief Description	Inhibitors
Purification	98.00%
Target Name	Carbonic Anhydrase inhibitor; Casein Kinase inhibitor; Hedgehog/Smoothened antagonist; PKA inhibitor; PKC inhibitor; Syk antagonist
Calculated MW	302.19
Formulation	C14H6O8
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
 SMILES: OC1=C(O)C2=C3C(=C1)C(=O)OC1=C3C(=CC(O)=C1O)C(=O)O2
 Pathways: GPCR/G Protein; Metabolism; Angiogenesis; Cytoskeletal Signaling; Stem Cells
 Receptor: Carbonic anhydrase; CK2; Smoothened; PKA; PKC; Syk
 Boling pt: 796.5
 Melting pt: 350
 Solubility: Ethanol:

References

1. Sekiguchi Y, et al. Bioorg Med Chem Lett. 2009 Jun 1;19(11):2920-3.

Note: This product is for in vitro research use only