

Decoquinatate

Catalog No: #S0471


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Decoquinatate
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	417.54
Formulation	C24H35NO5
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:ImmunologySMILES:CCCCCCCCCOc1c2c(c1)nc(=O)c(=O)c2OCCPathways:OthersReceptor:OthersBoiling pt:537.7CMelting pt:86-87CSolubility:DMSO: 0.01 mg/mL (0.02 mM)Appearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. The EFSA Journal, 2003, 17, 1-40.

Note: This product is for in vitro research use only