

Finasteride

Catalog No: #S0488


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Finasteride
Brief Description	Inhibitors
Purification	99.49%
Target Name	Reductase inhibitor
Calculated MW	372.55
Formulation	C23H36N2O2
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Endocrine

systemSMILES:C[C@]12CC[C@H]3[C@H]([C@@H]1CC[C@@H]2C(=O)NC(C)(C)C)CC[C@@H]4[C@@]3(C=CC(=O)N4)CPathways:Endocrinology /Hormones; EnzymeReceptor:5- α reductaseBoiling pt:Melting pt:253CSolubility:Ethanol: 75 mg/mL (201.32 mM); DMSO: 75 mg/mL (201.32 mM)Appearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Azzolina B, et al. J Steroid Biochem Mol Biol, 1997, 61(1-2), 55-64.

Note: This product is for in vitro research use only