

Tyrosine


Catalog No: #S0493

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Tyrosine
Brief Description	Inhibitors
Purification	99.00%
Target Name	Amino Acids and Derivatives inhibitor
Calculated MW	181.19
Formulation	C ₉ H ₁₁ NO ₃
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area:Endocrine systemSMILES:C1=CC(=CC=C1C[C@@H](C(=O)O)N)OPathways:OthersReceptor:Amino Acids and DerivativesBoiling pt:385.2CMelting pt:>300CSolubility:Water: 453 mg/L Appearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Dong S, et al. Blood. 2007 Jul 1;110(1):360-9.

Note: This product is for in vitro research use only