

Tetramisole hydrochloride


Catalog No: #S0494

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Tetramisole hydrochloride
Brief Description	Inhibitors
Purification	99.00%
Target Name	Dehydrogenase inhibitor
Calculated MW	240.75
Formulation	C ₁₁ H ₁₂ N ₂ SHCl
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Immunology
SMILES: Cl.C1CN2CC(N=C2S1)C1=CC=CC=C1
Pathways: Metabolism
Receptor: SDH
Boiling pt: 344.4°C
Melting pt: 266-267°C
Solubility: Water: 41 mg/mL (170.3 mM)
Appearance: White Crystal
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Eyre P, et al. J Pharm Pharmacol, 1970, 22(1), 26-36.

Note: This product is for in vitro research use only