

Creatinine

Catalog No: #S0496


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Creatinine
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	113.12
Formulation	C4H7N3O
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Metabolism systemSMILES:CN1CC(=O)N=C1NPathways:OthersReceptor:OthersBoiling pt:Melting pt:183CSolubility:Water: 23 mg/mL (203.32 mM)Appearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Allen PJ, et al. Neurosci Biobehav Rev, 2012, 36(5), 1442-1462.

Note: This product is for in vitro research use only