

Lauroseptol

Catalog No: #S0497


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Lauroseptol
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	364.46
Formulation	C19H42BrN
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OtherSMILES:[Br-].CCCCCCCCCCCCCCCC[N+](C)(C)CPathways:OthersReceptor:OthersBoiling pt:Melting pt:248-251CSolubility:Ethanol: 100 mg/mL (274.38 mM); Water: 100 mg/mL (274.38 mM); DMSO: 5 mg/mL (13.71 mM)Appearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Ito E, et al. Mol Pharmacol, 2009, 76(5), 969-983.

Note: This product is for in vitro research use only