

Oxybenzone


Catalog No: #S0499

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Oxybenzone
Brief Description	Inhibitors
Purification	99.92%
Target Name	Others
Calculated MW	228.24
Formulation	C14H12O3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: COC1=CC(O)=C(C(=O)C1=CC=CC=C1)
Pathways: Others
Receptor: Others
Boiling pt: 150-160°C
Melting pt: 62-64°C
Solubility: Water: 69 mg/L . Appearance: White Crystal
Remark: For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Gonzalez H, et al. Br J Dermatol. 2006 Feb;154(2):337-40.

Note: This product is for in vitro research use only