

Capitrol

Catalog No: #S0500


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Capitrol
Brief Description	Inhibitors
Purification	99.47%
Target Name	Others
Calculated MW	214.05
Formulation	C ₉ H ₅ Cl ₂ NO
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Immunology SMILES: OC1=C(Cl)C=C(Cl)C2=C1N=CC=C2 Pathways: Others Receptor: Others Boiling pt: Melting pt: 178-180°C Solubility: DMSO: 43 mg/mL (200.88 mM) Appearance: Gray Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Malaveille C, et al. IARC Sci Publ. 1991;(115):261-6.

Note: This product is for in vitro research use only