

Leucobasal

Catalog No: #S0501


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Leucobasal
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others inhibitor
Calculated MW	124.14
Formulation	C7H8O2
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Metabolism systemSMILES:COC1=CC=C(O)C=C1Pathways:OthersReceptor:radical polymerizationBoiling pt:243CMelting pt:56CSolubility:Ethanol: 25 mg/mL (201.38 mM); Water: 25 mg/mL (201.38 mM); DMSO: 25 mg/mL (201.38 mM)Appearance:Transparent Block CrystalRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Wang YH, et al. J AOAC Int. 2015 Jan-Feb;98(1):5-12.

Note: This product is for in vitro research use only