

Meglumine


Catalog No: #S0503

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Meglumine
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	195.21
Formulation	C7H17NO5
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OtherSMILES:CNC[C@@H]([C@H]([C@@H]([C@@H](CO)O)O)O)OPathways:OthersReceptor:OthersBoiling pt:490.4CMelting pt:129-131.5CSolubility:Water: 39 mg/mL (199.78 mM); DMSO: 39 mg/mL (199.78 mM)Appearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Navaei A, et al. Drug Dev Ind Pharm. 2014 Jun;40(6):701-10.

Note: This product is for in vitro research use only