

L-(-)- α -Methyldopa


Catalog No: #S0505

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	L-(-)- α -Methyldopa
Brief Description	Inhibitors
Purification	99.86%
Target Name	Dopamine Receptor inhibitor
Calculated MW	211.21
Formulation	C ₁₀ H ₁₃ NO ₄
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Cardiovascular System
 SMILES: C[C@](C)(CC1=CC(=C(C=C1)O)O)(C(=O)O)N
 Pathways: GPCR/G Protein; Neuroscience
 Receptor: DOPA decarboxylase
 Boiling pt: Melting pt: 300°C
 Solubility: Ethanol:

References

1. Head GA. Ann N Y Acad Sci. 1999 Jun 21;881:279-86.

Note: This product is for in vitro research use only