

4-Aminobutyric acid


Catalog No: #S0508

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	4-Aminobutyric acid
Brief Description	Inhibitors
Purification	99.00%
Target Name	GABA Receptor agonist
Calculated MW	103.12
Formulation	C4H9NO2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Nervous system
SMILES: NCCCC(O)=O
Pathways: Membrane transporter/Ion channel; Neuroscience
Receptor: GABA
Boiling pt: Melting pt: 195°C
Solubility: Soluble in DMSO
Appearance: White Crystal
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Erasso D, et al. Exp Neurol. 2014 Nov;261:646-53.

Note: This product is for in vitro research use only