

Dermofongin

Catalog No: #S0509


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Dermofongin
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	179.6
Formulation	C ₉ H ₆ ClNO
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Immunology SMILES: OC1=C2N=CC=CC2=C(Cl)C=C1 Pathways: Others Receptor: Others Boiling pt: Melting pt: 122-124°C Solubility: DMSO: 10 mM Appearance: Gray Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Latosiska JN, et al. J Mol Model. 2011 Jul;17(7):1781-800.

Note: This product is for in vitro research use only