

Iodoquinol

Catalog No: #S0513


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Iodoquinol
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	396.96
Formulation	C ₉ H ₅ I ₂ NO
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Immunology SMILES: OC1=C(I)C=C(I)C2=C1N=CC=C2 Pathways: Others Receptor: Others Boiling pt: 402°C Melting pt: >200°C Solubility: DMSO: 10 mM Appearance: Yellow Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Fisher AK, et al. J Toxicol Clin Toxicol. 1993;31(1):113-20.

Note: This product is for in vitro research use only