

Sodium dehydrocholate

Catalog No: #S0516


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Sodium dehydrocholate
Brief Description	Inhibitors
Purification	98.00%
Target Name	Lipid
Calculated MW	424.51
Formulation	C ₂₄ H ₃₃ NaO ₅
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Metabolism

systemSMILES: C[C@H](CCC(=O)[O-])[C@H]1CC[C@@H]2[C@@]1(C(=O)C[C@H]3[C@H]2C(=O)C[C@H]4[C@@]3(CCC(=O)C4)C)C.[Na+]Pathways: Metabolism Receptor: lipolysis Boiling pt: 581.5°C Melting pt: >300°C Solubility: DMSO: 10 mM Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Yousef IM, et al. Hepatology. 1987 May-Jun;7(3):535-42.

Note: This product is for in vitro research use only