

## Pancreatin

Catalog No: #S0517


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Pancreatin
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	N/A
Formulation	N/A
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:Metabolism systemSMILES:N/APathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Soluble in WaterAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Caruso C, et al. J Virol Methods. 2014 Nov;208:79-84.

Note: This product is for in vitro research use only