

Euquinine

Catalog No: #S0519


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Euquinine
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	396.48
Formulation	C ₂₃ H ₂₈ N ₂ O ₄
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Immunology SMILES: CCOC(=O)O[C@@H]([C@@H]1C[C@@H]2CCN1C[C@@H]2C=C)C3=C4C=C(C=CC4=NC=C3)OC Pathways: Others Receptor: Others Boiling pt: Melting pt: Solubility: DMSO: 10 mM Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

N/A

Note: This product is for in vitro research use only