

Diammonium Glycyrrhizinate

Catalog No: #S0522


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Diammonium Glycyrrhizinate
Brief Description	Inhibitors
Purification	97.00%
Target Name	Dehydrogenase inhibitor
Calculated MW	857
Formulation	C42H62O16N3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Metabolism

system SMILES: C[C@]12CC[C@](C[C@H]1C3=CC(=O)[C@H]4[C@]5(CC[C@H](C([C@H]5CC[C@]4([C@]3(CC2)C)C)(C)C)O[C@H]6[C@@H]([C@H]([C@@H]([C@H](O6)C(=O)O)O)O)[C@H]7[C@@H]([C@H]([C@@H]([C@H](O7)C(=O)O)O)O)C)(C)C(=O)O.N

Pathways: Metabolism
 Receptor: 11-beta-hydroxysteroid dehydrogenase
 Boiling pt: Melting pt: Solubility: Soluble in Water
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Jin J, et al. Food Chem Toxicol. 2014 Nov;73:95-104.

Note: This product is for in vitro research use only