

## 2,3,4-Trimethoxybenzaldehyde

Catalog No: #S0527


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	2,3,4-Trimethoxybenzaldehyde
Brief Description	Inhibitors
Purification	98%+
Target Name	Others
Calculated MW	196.2
Formulation	C10H12O4
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Other SMILES: COC1=C(OC)C(OC)=C(C=O)C=C1 Pathways: Others Receptor: Others Boiling pt: Melting pt: Solubility: Limited solubility Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

N/A

Note: This product is for in vitro research use only