

4-Fluoroanisole

Catalog No: #S0528


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	4-Fluoroanisole
Brief Description	Inhibitors
Purification	99%+
Target Name	Others
Calculated MW	126.13
Formulation	C7H7FO
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: COC1=CC=C(F)C=C1
Pathways: Others
Receptor: Others
Boiling pt: Melting pt: Solubility: Limited
solubility
Appearance: Liquid
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Laurence C, et al. J Phys Chem B. 2015 Feb 19;119(7):3174-84.

Note: This product is for in vitro research use only