

Dulcitol

Catalog No: #S0529


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|-------------------|---|
| Product Name | Dulcitol |
| Brief Description | Inhibitors |
| Purification | 98%+ |
| Target Name | Others |
| Calculated MW | 182.17 |
| Formulation | C6H14O6 |
| Storage | 3 years -20°C powder; 2 years -80°C in solvent; |

Images



Product Description

Research Area: Nervous system
SMILES: C([C@H]([C@@H]([C@@H]([C@H](CO)O)O)O)O)O
Pathways: Others
Receptor: Others
Boiling pt: Melting pt: Solubility: Soluble in Water
Appearance: White Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Jeon BS, et al. Enzyme Microb Technol. 2013 Aug 15;53(3):143-51.

Note: This product is for in vitro research use only