

Flavanone


Catalog No: #S0530

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Flavanone
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	224.26
Formulation	C15H12O2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: O=C1CC(OC2=CC=CC=C2)C1=CC=CC=C1
Pathways: Others
Receptor: Others
Boiling pt: Melting pt: Solubility: DMSO: 10 mM
Appearance: White Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Kean RJ, et al. Am J Clin Nutr. 2015 Mar; 101(3):506-14.

Note: This product is for in vitro research use only