

## Intermediate of Pregabalin


Catalog No: #S0531

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Intermediate of Pregabalin
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	187.24
Formulation	C <sub>9</sub> H <sub>17</sub> NO <sub>3</sub>
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Other SMILES: CC(C)C[C@H](CC(=O)N)CC(=O)O

Pathways: Others Receptor: Others Boiling pt: Melting pt: Solubility: Limited solubility Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

N/A

Note: This product is for in vitro research use only